





- · Reads both blood glucose and ketone test strips
- · Bluetooth connectivity
- Meets ISO 15197:20151
- · Downloads to Smartlog and Diasend/Glooko
- · Auto coding minimises errors · Easy to hold and grip
- · Stores 1,000 readings
- · Small sample size: 0.5 µl glucose/0.4 µl ketone
- · 5 second test time
- · Large, clear display
- · Strip ejection function

Blood Glucose Test Strips £9.95

SMARTLOG

1 Meter 2 Strips

Ketone Test Strips £9.95

Size 108 x 58 x 17 (mm) **Temperature** 5-45 °C

Relative humidity 10-90%

2015/02/10

Weight 72.6g (with batteries)

Item	Acceptance Criteria	Result		
Glucose Concentration <5.55mmol/L	95% of the measured glucose values falls within 0.83mmol/L	within <u>+</u> 0.83 mmol/L	within ±0.55 mmol/L	within <u>+</u> 0.28 mmol/L
		99.5%	96.6%	66.7%
Glucose	95% of the measured glucose values falls within <u>+</u> 15%	within <u>+</u> 15%	within <u>+</u> 10%	within <u>+</u> 1%
Concentration ≥5.55mmol/L		99.7%	95.7%	66.7%
99% of glucose values fall within zones A and B of the Consensus Error Grid		Within zones A and B 100%		

Approval number: HEA3868JUL21

References: 1. Data on file, i-sens ISO 15197:2015, system accuracy, May 2017



24/7 NON-AUTOMATED FREE SUPPORT 0800 881

Spirit Family of Companies Spirit House, Saffron Way Leicester, LE2 6UP

T: 0116 286 5000 F: 0116 271 9901

www.spirit-healthcare.co.uk